


WatchBP products positioning within the schema proposed by Dr. Pickering [6]



WatchBP Analyzer


WatchBP O3

Highly sophisticated software. Results are provided in PDF (see below) and Excel file format for analysis.

WatchBP home and O3

Easily transmits BP data from devices to PC for data collection, tabulation and analysis. User friendly interface.

Ambulatory Report of WatchBP O3



References

- World Health Organization. Guidelines for assessment and management of cardiovascular risk: 2007. www.ish-world.com/default.aspx?Guidelines
- O'Brien E, Asmar R, Belin L, et al. Practice guidelines of the European Society of Hypertension for clinic, ambulatory and self blood pressure measurement. *J Hypertens* 2005; 23: 697-701
- Parati G, Stergiou GS, Asmar R, et al. European Society of Hypertension guidelines for blood pressure monitoring at home: a summary report of the Second International Consensus Conference on Home Blood Pressure Monitoring. *J Hypertens* 2008;26:1505-26.
- The Japanese Society of Hypertension Guideline (JSH 2009). *Hypertension Research* (2009) 32, 11-23, 2009
- Mancia G, De Backer G, Dominiczak A, et al. 2007 Guidelines for the Management of Arterial Hypertension: The Task Force for the Management of Arterial Hypertension of the European Society of Hypertension (ESH) and of the European Society of Cardiology (ESC). *J Hypertens* 2007;25:1105-87.
- Pickering TG, Miller NH, Oggedegbe G, Krakoff LR, Artinian TN, Goff D, American Society of Hypertension; Preventive Cardiovascular Nurses Association. Call to action on use and reimbursement for home blood pressure monitoring: executive summary: a joint scientific statement from the American Heart Association, American Society of Hypertension, and Preventive Cardiovascular Nurses Association. *Hypertension* 2008; 52:1-9

For more information, please visit: www.watchbp.com

microlife

Be part of the great movement to prevent stroke and heart disease

WatchBP®

Be part of the great movement to prevent stroke and heart disease.

Let's work together for better health.



microlife

Microlife unique, guideline embedded, clinically validated BPM solution
WatchBP makes the great movement feasible & implementable.

Population screening

| Places | Features and Benefits |
|---|---|
| 1. Dr. office 2. Pharmacy 3. Work site 4. City hall ...etc. | Comprehensive solution with <ul style="list-style-type: none"> • Screen mode for identifying HTN (Hypertension) IAD (Inter Arm Difference) [1 -4, 6] AFIB (Atrial Fibrillation) • Routine mode: for easily implementing BP measurement • Auscultatory mode |



Diagnosis

| Places | Features and Benefits |
|------------|--|
| Dr. Office | <ul style="list-style-type: none"> • Smart 24 hrs ABPM with comprehensive software analyzer • Diagnostic mode for 7 days scheduling with readings tabulating in easy to read format [2 - 4, 6] |
| Dr. Office | <ul style="list-style-type: none"> • Diagnostic mode for 7 days scheduling with readings tabulating in easy to read format [2 - 4, 6] |



Treatment titration and patient education

| Places | Features and Benefits |
|------------|---|
| Dr. Office | <ul style="list-style-type: none"> • Selectable treatment target to fit patients with different conditions [4 - 6] • Display with color and gap to target for patient's education |



Treatment follow up & patient's home monitoring

| Places | Features and Benefits |
|----------------|--|
| Patient's home | <ul style="list-style-type: none"> • Diagnostic mode for 7 days scheduling with readings tabulating in easy to read format [2 - 4, 6] • Usual Mode for patients' home monitoring |



Early alert for high BP and Atrial Fibrillation

| Places | Features and Benefits |
|-------------------|--|
| Individual's home | <ul style="list-style-type: none"> • Simultaneously detecting high BP and AFIB conditions at home • Animated reminder for seeking medical consultation |

